




TUBERCULOSIS INFECTION T-CELL DETECTION SYSTEM


Simple operation, economical and convenient

Tuberculosis infection




Detection principle

After the human body is infected with Mycobacterium tuberculosis, sensitized T lymphocytes are produced and exist in the peripheral blood. When they are exposed to tuberculosis antigens (such as ESAT6, CFP10) again, they can be activated into effector T lymphocytes and release high-level cytokines. Among them, the most important one is interferon- γ (Interferon γ , IFN- γ). By detecting the concentration of interferon- γ , the tuberculosis infection situation can be judged. That is, the interferon- γ release assay, IGRAs.



Composition of tuberculosis infection T cell detection system



(1) Protein extraction reagent Stimulation and culture section



(2) Chemical Luminescence Platform Detection (MQ60, C800)

- Mono test CLIA(MQ60 series)
- Automatic CLIA (C800)

Detect the content of IFN - γ by chemiluminescence.



Stimulation and culture section
(protein extraction reagent)



Detection reagent
Quantitative detection * 3



N Tube: Blank Control Tube

Significance:

- 1) Eliminate the influence of γ interferon produced by the background.



T-tube: Test extraction tube, containing Mycobacterium tuberculosis-specific recombinant antigen (ESAT6).

Significance:

- 1) Specific antigens stimulate the sample to produce IFN- γ .




P - tube: Positive control tube, containing phytohemagglutinin (PHA) to induce lymphocytes to produce interferon.

Significance:

- 1) Evaluate whether the impact is caused by the low immunity of the sample itself (few lymphocytes).
- 2) Evaluate the impact caused by sample placement, incubation or operational techniques.

The process of γ -interferon in vitro release test



By comparing the IFN- γ release levels in the antigen-stimulated tube and the blank control tube, we can indirectly determine whether the subject has been infected with Mycobacterium tuberculosis.

Interpretation of γ -Interferon Release Assay Results

The indications of IGRA mainly include the diagnosis of latent tuberculosis infection and the auxiliary diagnosis of active tuberculosis.

Interpretation of IGRA test results:

- (1) A positive result supports the determination of the Mycobacterium tuberculosis infection status.
- (2) A negative result does not support the determination of the infection status.
- (3) An indeterminate result is one of the application limitations of this test.

Product Information of Hotgen Biotech

(1) Detection equipment



Analyzer Type	Fully automated Bench-top CLIA analyzer (MQ60 smart)
Sample type	Human serum, plasma, whole blood, urine, cerebrospinal fluid
Host size	300mm x 565mm x 610mm (lengthxwidthxheight)
Host weight	36.0kg
Languages	English
Mode	Batch mode
Channel	6 channels
Throughput	24 tests/hour
Storage	Result storage >1 million
Standby	24 hours operation standby
Data management	LIS/HIS connection

Detection Principle	AMPPD-Magnetic Particle Chemiluminescence
Throughput	150T/H
Sample Positions	30
Reagent Positions	12
First Report Time	17 min
Sample Types	Whole blood, Serum, plasma, urine
Sample Loading Needle	Steel needle
Sample Remove Needle	Steel needle
Dimensions (L x W x H)	520 x 760 x 650 mm



C800

(2) Detection reagent

Panel	Parameter	Qualification	Sample Type
IGRA	Protein extraction reagent	CE IVDR	WB
	Interferon- γ Chemiluminescence Immunoassay Kit (CLIA)	NMPA, CE	S/P

Simple operation
Economical and convenient
Be widely applied