

# Hotgen



## DENGUE NS1 Antigen Rapid Test (Colloidal Gold)

NS1 (non-structural protein 1) is a highly conserved glycoprotein that is essential for the viability of Dengue virus.

NS1 antigen is present with a high concentration in the blood of dengue patients during the early clinical phase of the disease, without being affected by previous infection history.

### NS1 Antigen Test

- *Excellent specificity*
- *High sensitivity*
- *Quick testing*
- *Simple operations*
- *Long shelf life*

 **Economy**

 **Rapid**

 **Precise**

 **Reliable**





## SPECIFICATION

Package	20T/kit
Sample Type	whole blood, serum or plasma
Test Time	15 min
Operation Temperature	10-30°C
Storage Temperature	2-30°C
Shelf Life (Unopened)	24 months

## PROCEDURES AND INTERPRETATION

